

iCool MediCube

Isothermic bag

for the transport of sensitive medication at an optimum temperature

▶ THE FREEDOM TO TRAVEL

Designed for the transport of sensitive medications at a temperature between **2 to 8°C for 36 hours**

▶ COOLING COMPARTMENT

keeps all your sensitive medication at a constant temperature

▶ STORAGE SPACE

For the transport of large amounts of medication over long distances. Cavity size of 9cm (width) x 18cm (length) x 14cm (height).

▶ RECOMMENDED USES

Insulin, growth hormones, adrenaline, monoclonal antibodies, vaccines, interferons

▶ MAXIMUM PROTECTION

Dimensions : 200 x 200 x 200 mm

Weight : 2 kg

