

iCool *Prestige*

Isothermic bag

for the transport of sensitive medication at an optimum temperature

▶ **THE FREEDOM TO TRAVEL**

Designed for the transport of sensitive medications at a temperature between 2 to 8°C for 24 hours

▶ **COOLING COMPARTMENT**

keeps all your sensitive medication at a constant temperature

▶ **STORAGE SPACE**

A special compartment reserved for your glucose (or other) meter, your syringes and needles

▶ **RECOMMENDED USES**

Insulin, growth hormones, adrenaline, monoclonal antibodies, vaccines, interferons

▶ **DISCREET AND ELEGANT**

Dimensions : 190 x 170 x 60 mm

Weight : 780 g

