

EasyBag

Single

Isothermic bag

for the transport of insulin or adrenaline at an optimal temperature

▶ **THE FREEDOM TO TRAVEL**

Designed for the transport of one pre-filled insulin or adrenaline pen at an optimum temperature between 15 to 25 ° C for 3 days

▶ **COOL WITHOUT ELECTRICITY**

The EasyBag Single becomes active in less than 40 seconds when the inner pouch is soaked in plain tap water. As long as the inner pouch is wet, your medication is kept at the right temperature

▶ **PATENTED TECHNOLOGY BY BASF IN GERMANY**

▶ **CAN BE REUSED HUNDREDS OF TIMES**

QUICK AND EASY TO USE

▶ **IDEAL FOR ACTIVE TRAVELLERS**

▶ **RECOMMENDED USES**

▶ Pre-filled pens of insulin or adrenaline

DISCREET AND ELEGANT

Dimensions : 210 x 55 x 30 mm

Poids : 70 g

