

iCool MediCube ULTRA

Isothermal bag

for the transport of sensitive medication
between the temperatures of 2 and 8°C



iCool MediCube ULTRA features

THE FREEDOM TO TRAVEL

Designed for the transport of sensitive medications at a temperature between 2 to 8°C for 36 hours

COOLING COMPARTMENT

keeps all your sensitive medication at a constant temperature.
Includes 3 x 900g polymer gel packs

STORAGE SPACE

For the transport of large amounts of medication over long distances.
Useable space of 21 x 14 x 12.5cm (W x L x H)

RECOMMENDED USES

Insulin, growth hormones, adrenaline, monoclonal antibodies, vaccines, interferons, Anti-TNF medications


Specially formulated polymer that stays frozen twice as long as a traditional ice packs

Holds up to 40 injection pens or syringes



36 hours
of protection
between 2 and 8°C

The iCool MediCube ULTRA keeps all your medication at a temperature between 2 and 8°C for 36 hours. It is the perfect travel companion for longer trips.



- Weight: 3kg
- Dimensions: 260 x 260 x 160mm
- Cavity Size: 210 x 140 x 125
- Operating Environment:
Temp: 5 °C ~ 40 °C, humidity: ≤90% RH
- Conditions for transport and storage:
Temp: -20 °C ~ 60 °C, humidity <93% RH

Discover all the MedActiv solutions
for the safe transport of medication
www.medactiv.com.au

Imagined and developed in France by
HealthWorks France
75008, Paris, France
contact@medactiv.com



medactiv
transforming lives